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| Product                                    | Raw Material                        | Standard           |
|--------------------------------------------|-------------------------------------|--------------------|
| Healing Abutment                           | Titanium-6Al-4V ELI                 | ASTM F136          |
| Cover Screw                                | Titanium-6Al-4V ELI                 | ASTM F136          |
| Combi Abutment                             | Pure titanium grade 4               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Dual Abutment                              | Pure titanium grade 4               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Custom Abutment                            | Pure titanium grade 4               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Dual Milling Abutment                      | Pure titanium grade 4               | ASTM F67           |
|                                            | Pure titanium grade 2               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Angled Abutment                            | Pure titanium grade 4               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Ti-Temporary Abutment                      | Pure titanium grade 4               | ASTM F67           |
|                                            | Pure titanium grade 2               | ASTM F67           |
| Screw Abutment                             | Pure titanium grade 4               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Angled Screw Abutment                      | Pure titanium grade 4               | ASTM F67           |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |
| Ti-Cylinder                                | Titanium-6Al-4V ELI                 | ASTM F136          |
| Ti-Temporary Cylinder                      | Pure titanium grade 2               | ASTM F67           |
| Direct casting Abutment,<br>Gold Cylinder  | Gold Alloy                          | N/A                |
|                                            | Polyacetal (POM)                    | CAS No. 66455-31-0 |
| Metal Casting Abutments,<br>Metal Cylinder | Cobalt-28Chromium-6Molybdenum Alloy | ASTM F799          |
|                                            | Polyacetal (POM)                    | CAS No. 66455-31-0 |
| Scan Abutment                              | Polyetheretherketone (PEEK)         | ASTM F2026         |
|                                            | Titanium-6Al-4V ELI                 | ASTM F136          |

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