Dentium
For Dentists By Dentists


Contents

Lab Products

- 04 rainbow™ Multi-Layer Block
- 07 rainbow™ Indication
- 08 rainbow™ Shine-T Block
- 12 rainbow™ High Shine Block
- 15 rainbow™ Shade Block
- 17 rainbow™ Trans Block
- 19 rainbow™ Block
- 22 rainbow™ LS Block
- 24 rainbow™ Shine Liquid
- 31 rainbow™ Coloring Guide
- 33 rainbow™ Paste Stain
- 40 rainbow™ Paste Stain Guide
- 49 rainbow™ Porcelain
- 55 rainbow™ Porcelain Guide

Digital Instrument

- 57 rainbow™ Mill-Metal
- 59 rainbow™ Mill - Zr 2nd
- 62 rainbow™ Mill Bur List

rainbow™
Block System

By Dr. Ha-Young Kim


rainbow™ Multi-Layer Block

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
rainbow™ Multi-layer Block Indication

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




A0

A2

rainbow™ Shine-T Block Clinical Case 1

Brushing

Final prosthesis

rainbow™ Shine-T Block Clinical Case 2

Brushing

Final prosthesis

rainbow™ Shine-T Block Clinical Case 3

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


rainbow™ High Shine Block Clinical Case 1 - Full Zirconia

CAD design

Brushing


Sintering

Final prosthesis


rainbow™ High Shine Block Clinical Case 2 - Full Zirconia

CAD design

Brushing


Sintering

Final prosthesis

rainbow™ High Shine Block Clinical Case 3 - Full Zirconia

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

A0.5

A2


rainbow™ Shade Block Clinical Case 1 - Build-up

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




rainbow™ LS Block Clinical Case

By Dr. Hyun-Sik Park _ HEIM Dental Clinic

Crystallization를 위한 Firing parameter

Stand-by temperature B[°C]	Closing time S[min]	Heating rate t1[°C/min]	Firing temperature T1[°C]	Holding time H1[min]	Heating rate t2[°C/min]
403	6:00	90	820	0:10	30

Firing temperature T2[°C]	Holding time H2[min]	Vacuum1 1 ₁ [°C]/1 ₂ [°C]	Vacuum2 2 ₁ [°C]/2 ₂ [°C]	Long-term cooling L[°C]	Cooling rate tl[°C/min]
840	7:00	550/820	820/840	700	0

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* Shine Liquid Pen

Type	Volume	Art. No.				
Bottle	20 ml	PBA1	PBA2	PBA3	PBA3.5	PBA4
		PBB1	PBB2	PBB3		PBB4
		PBC1	PBC2	PBC3		PBC4
			PBD2	PBD3		PBD4
	100 ml	PBA1	PBA2	PBA3	PBA3.5	PBA4
		PBB1	PBB2	PBB3		PBB4
		PBC1	PBC2	PBC3		PBC4
			PBD2	PBD3		PBD4

Type	Volume	Art. No.		
Pen	—	PEN-TV (Trans Violet)	PEN-LG (Light Grey)	PEN-DG (Dark Grey)
Bottle	20 ml	TV (Trans Violet)	LG (Light Grey)	DG (Dark Grey)
	100 ml	TV (Trans Violet)	LG (Light Grey)	DG (Dark Grey)


Type	Volume	Art. No.
Bottle	20 ml	P-Red (지시색소)
	100 ml	Thinner (희석액)

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



Sintering Schedule

Shine-T Block/ High Shine Block/ Shade Block/ Trans Block

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







Use example of rainbow™ Paste Stain

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Use example of rainbow™ Paste Stain

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rainbow™ Veneering & Staining Material

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



Physical & Chemical Properties of rainbow™ Porcelain

The results of physical and chemical properties of rainbow™ Porcelain are measured in accordance with ISO and KFDA standard, and suitable to all of these standards.

Properties	Unit	Results	ISO standard
CTE ¹⁾ (25~500°C)	10 ⁻⁶ /°C	9.3	± 0.5
Glass Transition Temperature	°C	570	± 20
Chemical Solubility	μg/cm ²	1.5	≤ 100
Radioactivity	Bq/g	0.0003	≤ 1.0
Average Particle Size	μm	30±3	-

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Shade	White	Trans
Fluor Effect	FLUOR-W	FLUOR-T

Shade	Body	Trans
Add-on	ADD-B	ADD-T

Glaze & Stain	GLAZE	ST-Yellow	ST-Ocher	ST-Orange	ST-Brown	ST-yBrown	ST-rBrown	ST-dBrown	ST-Pink
---------------	-------	-----------	----------	-----------	----------	-----------	-----------	-----------	---------

ST-Lilac	ST-sBlue	ST-cBlue	ST-Black	ST-Grey	ST-White	ST-A	ST-B	ST-C
----------	----------	----------	----------	---------	----------	------	------	------

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


Special Liquids of rainbow™ Porcelain

Liquid	Unit	Description
Liner Liquid	10ml, 100ml	Liquid for painting Liner porcelain uniformly on zirconia coping
Build-Up Liquid	10ml, 100ml	Liquid for Body build-up
Margin Liquid	10ml, 100ml	Liquid for making sharp edge of Margin
Glazing Liquid	10ml, 100ml	Liquid for painting Glaze and Stain uniformly

Kit of rainbow™ Porcelain

Professional Kit

Type	Color	Unit	Contents
Liner	7	10g	LNB, LNT, LN1, LN2, LN3, LN4, LN5
Opacious Dentin	17	10g	ODA1, ODA2, ODA3, ODA3.5, ODA4, ODB1, ODB2, ODB3, ODB4, ODC1, ODC2, ODC3, ODC4, ODD2, ODD3, ODD4, ODW
Dentin	20	10g	DA1, DA2, DA3, DA3.5, DA4, DB1, DB2, DB3, DB4, DC1, DC2, DC3, DC4, DD2, DD3, DD4, DA00, DA0, DB00, DB0
Enamel	3	10g	EN1, EN2, EN3
Translucent	8	10g	TNT, TYL, TOR, TBR, TPK, TBL, TGR, TML
Opal & Fluor Effect	6	10g	OTR, OE1, OE2, OE3, FLUOR-W, FLUOR-T
Dentin Modifier	7	10g	DM-YL, DM-OR, DM-BR, DM-PK, DM-BL, DM-GR, DM-WT
Cervical	4	10g	CVA, CVB, CVC, CVD
Margin	5	10g	MG1, MG2, MG3, MG4, MG5
Gum	3	10g	GLP, GNP, GDP
Add-on	2	10g	ADD-B, ADD-T
Liquid	4	10ml	Liner Liquid (1), Build-Up Liquid (2), Margin Liquid (1)

Advanced Kit

Type	Color	Unit	Contents
Liner	4	10g	LNB, LN1, LN2, LN3
Opacious Dentin	4	10g	ODA2, ODA3, ODB2, ODB3
Dentin	4	10g	DA2, DA3, DB2, DB3
Enamel	4	10g	EN2, EN3, OE2, OE3
Translucent	4	10g	TNT, TOR, TBL, TML
Cervical	2	10g	CVA, CVB
Margin	2	10g	MG2, MG3
Glaze & Stain	1	3g	Glaze
Liquid	4	10ml	Liner Liquid (1), Build-Up Liquid (1), Margin Liquid (1), Glazing Liquid (1)

Intro Kit

Type	Color	Unit	Contents
Liner	2	10g	LNB, LN2
Opacious Dentin	1	10g	ODA2
Dentin	1	10g	DA2
Enamel	2	10g	EN2, OE2
Translucent	2	10g	TNT, TOR
Glaze & Stain	1	3g	GLAZE
Liquid	2	10ml	Liner Liquid, Build-Up Liquid


Stain Kit

Type	Color	Unit	Contents
Glaze & Stain	18	3g	GLAZE, ST-Yellow, ST-Ocher, ST-Orange, ST-Brown, ST-yBrown, ST-rBrown, ST-dBrown, ST-Pink, ST-Lilac, ST-sBlue, ST-cBlue, ST-Black, ST-Grey, ST-White, ST-A, ST-B, ST-C
Liquid	1	10ml	Glazing Liquid

Final

rainbow™
Porcelain Guide

Labial

Lingual

Lateral

Use example of rainbow™ Porcelain

Pre-Liner

Liner : LN3, LNT

Applying as thin as possible
but avoiding flow into Coping inside

Fluor-T

Applying thin layer on cervical for expressing
Natural Tooth's Fluorescence effect

Dentin : DA3

Giving natural tooth's structure & color

Dentin : DA2

Giving natural tooth's structure & color

Translucent : TOR

Partially applying on Light's reflected area

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
TOR : Partially applying on Light's reflected area
TBL : Giving Transparent Effect
OTR : Providing Opal Effect

Opal Effect : OE1, Enamel : EN2

OE1 : Improving Opal Effect & Brightness
EN2 : Expressing Translucent Effect on Incisal

Enamel : EN3

EN3 : Expressing Translucent Effect on Incisal

Dentin : DA2, DA3

EN2, EN3 : Expressing Translucent Effect on Incisal

Opal Effect : OE1

OE1 : Improving Opal Effect & Brightness

Final

Partially applying on Light's reflected area

Use example of rainbow™ Porcelain

After Firing

Contouring

After Glazing

rainbow™ Digital Instrument

Scan QR Code to
You can see Video


Optimized for implant products from Dentium


Dentium Ti Blank

Hex

Non-Hex

Cross section of Dentium Customized titanium Abutment & SuperLine FX4510 SW

Optimized for Glass ceramic, Zirconia, PMMA

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
Material	Premill Abutment, PMMA/Acryl, Glass ceramic, Hybrid Ceramic
File Load	LAN / Wi-fi
Control axis	4 Axis
Spindle	60,000rpm / Max 1200W
Motion drives	All servo motors
Tool Posts	8 Posts (Automatic tool change)
Ti Blank holder	10
Milling type	Wet processing
Dimensions	51 x 71 x 64 cm
Weight	150kg

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
Material	Zirconia, PMMA, Wax
Control axis	Simultaneous 5-axis
Spindle	60,000 rpm, 500W
CAM S/W	hyperDENT
Tool posts	7 Posts (Shank Ø4)
Milling type	Dry processing
Dimensions (W x D x H)	35 x 54 x 64 cm
Weight	90 kg

rainbow™ Mill Bur List

rainbow™ Mill - Zr 2nd

Design	Art. No.	Stopper Color	Shank diameter	Diameter	Material
	BE4D2-D	White	Ø4	Ø 2	ZIRCONIA
	BE4D1-D			Ø 1	
	BE4D0.5-D			Ø 0.5	
	BE3D2-N	Blue	Ø 4	Ø 2	PMMA / WAX
	BE3D1-N			Ø 1	

rainbow™ Mill - Metal

Design	Art. No.	Type	Shank diameter	Diameter	Material
	BE6D3-T	ENDMILL	Ø 6	Ø 3	TITANIUM
	BE6D2-T	ENDMILL	Ø 6	Ø 2	
	BE6D1.5-T	ENDMILL	Ø 6	Ø 1.5	
	BE6D1.0-T (옵션)	ENDMILL	Ø 6	Ø 1.0	
	BE6D2-N	ENDMILL	Ø 6	Ø 2	PMMA / WAX
	BE6D1-N	ENDMILL	Ø 6	Ø 1	
	GT6D2-ED	GRINDING	Ø 6	Ø 2	GLASS CERAMIC / COMPOSITE RESIN
	GT6D1-ED	GRINDING	Ø 6	Ø 1	
	GT6D0.6-ED	GRINDING	Ø 6	Ø 0.6	