

Guided Surgery Process

Dentium Digital Center



DC DIGITAL CENTER
www.dentiumdc.co.kr

DENTIUM GUIDE SYSTEM



Digital Guide

Simple & Accurate
Digital Guided Surgery

DENTIUM GUIDE SYSTEM

Polymer Guide

Thermoplastic Surgical Template

Implant Guide

Surgical guide utilizing silicon spacer
and unique parallel pin

Dentium
For Dentists By Dentists

HEAD OFFICE
서울시 강남구 테헤란로87길 21, 10층(삼성동, 동성빌딩) Tel 02-555-3750 Fax 02-501-9560

Specifications are subject to change without any notice.
Some products listed in this catalog may not be available in the market due to pending approval.

DGSL-1703 [Rev.1]



Dentium
For Dentists By Dentists

NEW
PRODUCT



Digital Guide

Simple & accurate digital guided surgery



Simple Kit

- Accurate Depth and Positions can be obtained in a narrow space
- Application of Ø3.0 sleeve hole, minimization of planning constraints due to adjacent tooth interference




Full Kit

- Progression from drilling to implant placement.
- Step drill shape securing guide section, ensuring convenience and accuracy of procedure.



Polymer Guide

Thermoplastic surgical template for dental implant placement




Construction planning and template production in plaster model



Single Standard




Cantilever Multi-Ready



Implant Guide

Surgical guide utilizing silicon spacer and unique parallel pin



Guide Drill + Silicon Spacer



Ellipse Parallel Pins



Tripod Parallel Pins