

bright *X-ray*

PORTABLE & SUITABLE



Dentium

bright X-ray



Patient Type

Adult, Child

Tooth Type

Anterior Tooth, Canine Tooth,
Molar Tooth and Premolar Tooth

Irradiation Time

Time and Angle can be controlled

Product Information



Size : 304.3 X 160 X 262.5 (mm)

Weight : 1.6kg (Main body)

Battery Specifications : 11.1V 1,500mAh

※ Available for taking picture of 200 X-ray
photographs with 2 hours' charging

Light weight Clear Image



Hand Strap



Wrist Rest



Post-scatter
Shielding Glass



Stand

Accessories

Dentium