

Precise with Guide Pin

Pin Guide *Simple*




Mapping
Alveolar Bone


Optimized
for Flapless Surgery

Excellent
Performance
in Narrow Ridge Case

Dentium



Pin Guide *Full*



Pin Guide System


Specification

XPGSK

(unit: mm)

Blade Marker

Use: Check drilling position



XPBM


Tip

Handle

SIZE 1:3

Guide Pin




XSGP0015




* Disposable

How to Use

Bend by **hand** according to the 'Bending Guide'




Procedure of Use




Clinical Case

I Case 1. Overcoming Narrow Ridge

Pre-OP



Post-OP




I Case 2. Simple Posterior Implantation with Flapless


Pre-OP



Post-OP




Pre-OP




Pin Insertion



Drilling




Implant Placement



Post-OP



Precise with Guide Pin



Dentium